



**BIOANTIBODY**

# Dengue NS1 Antigen Rapid Test Kit

(Lateral chromatography)

## About Dengue

Dengue viruses are spread to people through the bite of an infected *Aedes* species (*Ae. aegypti* or *Ae. albopictus*) mosquito. Almost half of the world's population, about 4 billion people, live in areas with a risk of dengue. Dengue is often a leading cause of illness in areas with risk.


*A blood test is the only way to confirm the diagnosis.*



## Intended Use

Dengue NS1 Antigen Rapid Test Kit (Lateral chromatography) is designed for the early detection of dengue virus NS1 antigen in human serum, plasma, whole blood or fingertip whole blood.

Immunochromatographic and double-antibody sandwich method is applied to detect Dengue NS1. This test is for professional use only.



## Key Advantages



### Stable & sufficient supply chain

self-supplied key raw materials



### Easy to operate

no instruments needed, visually read results



### Rapid to get results


just 15 mins



### Economic to store and transport

2~30°C, don't require any cold chain facilities

## Operation Flow



Note: do NOT read the results after 20 mins!

## Result Interpretation

Positive	Negative	Invalid
<p>Two lines (both T and C line) appear</p>	<p>Only the C line appears and the detection T line does not appear</p>	<p>No line appears in C line no matter T line appears or not</p>

## Order Information

Product Name	Cat. No	Size	Specimen	Shelf Life	Trans. & Sto. Temp.
Dengue NS1 Antigen Rapid Test Kit (Lateral chromatography)	B010C-01	1test/kit	S/P/WB	24 Months	2-30°C / 36-86°F
	B010C-25	25 tests/kit			





## BIOTECH FOR A BETTER LIFE



**BIOANTIBODY**

### Bioantibody Biotechnology Co., Ltd.

Website: [www.bioantibody.net](http://www.bioantibody.net)

Tel: +(86)025-58501357

E-mail: [Sales@bkbio.com.cn](mailto:Sales@bkbio.com.cn)

Address: Room 903/905, Building C6, 9 Weidi Road, Nanjing, Jiangsu, China

